

Service Report: 1016**Date:** August 25, 2026**Hospital / Client:** Joint Medical Company**Device:** Ems Radial Shockwave (Dolarclast)**Serial Number:** LQ3220**Installation Date:** February 01, 2026**Department:** Physiotherapy Dept**Issue Reported:**

End user reported weak shock output from the EMS Radial Shock Wave device,
The maintenance indicator light shown on the device.

Findings Upon Inspection:

Upon inspection, it was found that each shockwave kit reaches its usable life and produces weak shocks once the count limit is reached,
Device installed on 1-2-2026. Total shock count from installation date to present: 9,131,371 shocks.
Kit 3: reached 1,543,126 shocks, replaced on 25-8-2026,
Kit 2: reached end of life, replaced on 19-7-2026,
Kit 1: reached end of life, replaced on 02-3-2026.

Action Taken:

Replaced the old kit with a new one and performed a test run after each replacement to confirm normal shock output.

Status:

- Resolved

Recommendation

Track the shock count per kit regularly and replace proactively once it approaches the manufacturer's recommended limit, to avoid weak-shock complaints from the end user.

Customer Sign	Hospital B.M.E.	Orthoserv Biomedical
		Mahmoud Ghasan

© 2026 OrthoServ

ص.ب. 286869 الرياض 11323 المملكة العربية السعودية. عضوية رقم 34907 هاتف 014806641 فاكس 012419801
P.O.Box 286869 Riyadh 11323 Kingdom of Saudi Arabia Membership No 34907 Tel 014806641 Fax: 012419801

CR: 1010251830 ORTHOPEDIC SERVICE EST